

# Power Port<sup>\*</sup>

IMPLANTABLE PORT

Feel and See the **NEW**  
Standard of Care<sup>\*</sup>



**BAIRD**  
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.<sup>\*</sup>*



A portion of the proceeds of the sale of this device will be used to support breast cancer initiatives.



## It's About the Patient...

Upon receiving a **PowerPort\*** Implanted port, patients are provided a **Discharge Packet** to identify and educate patients which includes:



\* Standard for all PowerPort\* Implanted ports



## PowerLoc\* Safety Infusion Set:

The **PowerLoc\*** Safety Infusion Set is the ideal device for accessing the PowerPort\* device for power injection of contrast media.

### ○ Power Injectable:

- Withstands 300 psi, at Flow Rates up to 5 ml/s
- Small Footprint
- Easy to Dress
- Ultra Low-Profile



Feel and See the **NEW**  
Standard of Care\*

## Features & Benefits Unique to **PowerPort\* duo M.R.I.\*** Device:

**PowerPort\* duo M.R.I.\*** Dual Lumen Implanted Port is ideal for delivery of supportive therapies in conjunction with intravenous chemotherapy.

Intermediate in size, the PowerPort\* duo M.R.I.\* device requires only a small pocket to facilitate easy placement while ensuring greater patient comfort.

### ○ **Enables Reliable Venous Access**


- Reduction of Repeated Needlesticks
- Easy, Flexible Placement and Access
- Power Injectable for CECT Scans
- Delivery of I.V. Therapies

### ○ **Enables Superior Imaging**

- Power Injection and CECT Scan  
Compatible when used with PowerLoc\*  
Safety Infusion Set
- Withstands 5 ml/s flow rate @ 300 psi  
Pressure Limit Setting
- Power Injected CECT Scans Produce  
Superior Images to Help Better-Manage  
Patient Care

### ○ **Easily Identifiable**

- Soft Palpation Points on Septum
- Radiopaque Identifier
- Unique Patient Discharge Packet



Flipped Port



X-Ray Correct Identification

Unique Radiopaque Symbol for Power Identification and Flipped Port Detection.



**PowerPort\* duo M.R.I.\* Implanted Port - Radio Translucent**

Product Code	Kit Type	Code Description
1829500	Intermediate Kit	PowerPort* duo M.R.I.* Port w/ 9.5 Fr ChronoFlex* Catheter <b>with</b> Suture Hole Plugs
1829570	MicroIntroducer Kit	PowerPort* duo M.R.I.* Port w/ 9.5 Fr ChronoFlex* Catheter <b>with</b> Suture Hole Plugs



**PowerLoc\* Safety Infusion Set - with Y-Site**

Product Code	Code Description with Y-Site	Max. Recommended Flow Rate
0671934	PowerLoc* Safety Infusion Set, 19 Ga x .75"	5 ml/s
0671910	PowerLoc* Safety Infusion Set, 19 Ga x 1"	5 ml/s
0671915	PowerLoc* Safety Infusion Set, 19 Ga x 1.5"	5 ml/s
0672034	PowerLoc* Safety Infusion Set, 20 Ga x .75"	5 ml/s
0672010	PowerLoc* Safety Infusion Set, 20 Ga x 1"	5 ml/s
0672015	PowerLoc* Safety Infusion Set, 20 Ga x 1.5"	5 ml/s
0672234	PowerLoc* Safety Infusion Set, 22 Ga x .75"	2 ml/s
0672210	PowerLoc* Safety Infusion Set, 22 Ga x 1"	2 ml/s
0672215	PowerLoc* Safety Infusion Set, 22 Ga x 1.5"	2 ml/s

**PowerLoc\* Safety Infusion Set - without Y-Site**

Product Code	Code Description without Y-Site	Max. Recommended Flow Rate
0651934	PowerLoc* Safety Infusion Set, 19 Ga x .75"	5 ml/s
0651910	PowerLoc* Safety Infusion Set, 19 Ga x 1"	5 ml/s
0651915	PowerLoc* Safety Infusion Set, 19 Ga x 1.5"	5 ml/s
0652034	PowerLoc* Safety Infusion Set, 20 Ga x .75"	5 ml/s
0652010	PowerLoc* Safety Infusion Set, 20 Ga x 1"	5 ml/s
0652015	PowerLoc* Safety Infusion Set, 20 Ga x 1.5"	5 ml/s
0652234	PowerLoc* Safety Infusion Set, 22 Ga x .75"	2 ml/s
0652210	PowerLoc* Safety Infusion Set, 22 Ga x 1"	2 ml/s
0652215	PowerLoc* Safety Infusion Set, 22 Ga x 1.5"	2 ml/s

**Minimum Gravity Flow Rate for PowerLoc\* Safety Infusion Set**

19 Ga.	1,680 ml/hr
20 Ga.	960 ml/hr

**Maximum Priming Volume for PowerLoc\* Safety Infusion Set**

19 Ga., 20 Ga. and 22 Ga.	<b>with</b> Y-Site - 0.4 ml
19 Ga., 20 Ga. and 22 Ga.	<b>without</b> Y-Site - 0.3 ml



ACCESS SYSTEMS

**Bard Access Systems**

605 North 5600 West  
Salt Lake City, Utah 84116  
801-595-0700  
www.bardaccess.com  
www.powerportadvantage.com

Customer Service: 800-545-0890  
Clinical Information: 800-443-3385

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

\*Bard, PowerPort, PowerLoc, M.R.I., "Feel and See the New Standard of Care.", "Advancing the Delivery of Health Care." and the color Purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. ChronoFlex is a registered trademark of Cardiotech International, Inc. or an affiliate. Taxol is a registered trademark of Bristol-Myers Squibb Company. or an affiliate. Taxotere is a registered trademark of Aventis Pharma S. A. or an affiliate.

© Copyright 2009 C. R. Bard, Inc. All rights reserved.

MC-0408-00 0901R