

Power PICC^{*}
THE UNIVERSAL PICC

The Power of Purple^{*}



BAIRD
ACCESS SYSTEMS

Advancing the Delivery of Health Care.^{}*

The Power of Purple*

Features and Advantages:

- **Maximum injection rates –**
allow injection of contrast media for CECT scans at a maximum rate of 5 ml/sec
- **Easy identification –**
catheters are clearly labeled to differentiate them from all non-indicated power injection PICCs
- **Easy-to-read ID tags –**
indicate maximum injection rates and prompt users to check catheter patency
- **Reverse taper design –**
gently plugs the insertion site and provides superior kink resistance
- **CVP monitoring –**
all PowerPICC* catheters are indicated for CVP monitoring
- **Multiple wire configurations –**
each wire type is available in 70 cm and 135 cm wire lengths
- **Patient satisfaction –**
improves patient comfort due to decreased needle sticks


Multiple Configurations:

| Single-Lumen | Dual-Lumen | Triple-Lumen |
|--------------|------------|--------------|
| 4 Fr. | 5 Fr. | 5 Fr. |
| 5 Fr. | 6 Fr. | 6 Fr. |

Premium Guidewire Options:

| |
|---|
| Nitinol Guidewire - 70 cm & 135 cm |
| Stainless Steel Guidewire* - 70 cm & 135 cm |


*Not available in 5 Fr. Triple Lumen



Both the nitinol and stainless steel guidewires are available in 70 cm and 135 cm lengths



PowerPICC* catheters combine the efficacy of PICC access with the efficiency of power injection for CECT scans





Ordering Information

Product Specifications:

| Type of Catheter | Product Code | Guidewire | Gravity Flow (ml/hr) | Priming Volume (ml) | Gauge |
|------------------------------|--------------|------------------------|----------------------|---------------------|-------|
| 4 Fr. SL PowerPICC* Catheter | 3174135 | 135 cm – standard wire | 1,272 | 0.67 | 18 |
| | 3174155 | 70 cm – standard wire | | | |
| | 3174335 | 135 cm – nitinol wire | | | |
| | 3174355 | 70 cm – nitinol wire | | | |
| 5 Fr. SL PowerPICC* Catheter | 3175135 | 135 cm – standard wire | 1,185 | 0.66 | 18 |
| | 3175155 | 70 cm – standard wire | | | |
| | 3175335 | 135 cm – nitinol wire | | | |
| | 3175355 | 70 cm – nitinol wire | | | |
| 5 Fr. DL PowerPICC* Catheter | 3275135 | 135 cm – standard wire | 578/578 | 0.57/0.57 | 18/18 |
| | 3275155 | 70 cm – standard wire | | | |
| | 3275335 | 135 cm – nitinol wire | | | |
| | 3275355 | 70 cm – nitinol wire | | | |
| 6 Fr. DL PowerPICC* Catheter | 3276135 | 135 cm – standard wire | 753/753 | 0.62/0.62 | 18/18 |
| | 3276155 | 70 cm – standard wire | | | |
| | 3276335 | 135 cm – nitinol wire | | | |
| | 3276355 | 70 cm – nitinol wire | | | |
| 5 Fr. TL PowerPICC* Catheter | 3385335 | 135 cm – nitinol wire | 982/131 | 0.76/0.43 | 18/20 |
| | 3385355 | 70 cm – nitinol wire | | | |
| 6 Fr. TL PowerPICC* Catheter | 3386135 | 135 cm – standard wire | 1,163/275 | 0.76/0.47 | 17/19 |
| | 3386155 | 70 cm – standard wire | | | |
| | 3386335 | 135 cm – nitinol wire | | | |
| | 3386355 | 70 cm – nitinol wire | | | |

PowerPICC* Catheter Tray Components:

| Description |
|---|
| PICC |
| Microintroducer, 10 cm |
| Scalpel |
| Introducer needle, 21 Ga. |
| End cap (one per lumen) |
| Syringe, 12 ml luerlock, latex free |
| StatLock* catheter stabilization device |
| Measuring tape |
| Stainless steel guidewire/Nitinol guidewire |

The Power of Purple*

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Learn More - For more information on the PowerPICC* Catheter, please contact your Bard Access Systems representative or visit www.powerpicc.com. Call toll-free at 1-800-545-0890, or visit us on the web at www.bardaccess.com.



ACCESS SYSTEMS

Bard Access Systems, Inc.

Salt Lake City, Utah 84116 USA
801-595-0700
www.bardaccess.com

Customer Service: 800-545-0890
Clinical Information: 800-443-3385